


# Never clean a commode bucket again – for less than the price of a Sunday newspaper.


Medline’s disposable commode liners make commode cleanup easy and hygienic. Each liner includes an absorbent pad made from a polymer gelling agent that solidifies and neutralizes waste for convenient disposal. Because they are solidified, liquids will not move around in the commode, making transport easier. All you need to do is toss the liner in the garbage after use.

Disposable commode liners also eliminate contact between staff and human waste, to help reduce contamination. The average person will use two liners each day, at a suggest retail price of \$1.50 each.



Item	Description	Pkg
<b>Disposable Commode Liner</b>		
MDS89664LINER	Commode Liner	12 liners/box, 6 boxes/cs

## How to use Disposable Commode Liners



1. Place the commode liner in the commode. Make sure that the liner is folded over the edge of the bucket.
2. As liquid hits the absorbent core, the polymer inside the core immediately begins converting the liquid into a solid. The core neutralizes odor as the liquid solidifies.
3. Simply toss the liner into the garbage after use.