

# Guardian Easy Care® Commodes



Model G30213-4F

## 3 IN 1 STEEL COMMODE

---

- Folding frame design assembles in seconds without tools and reduces storage space by 35% while making delivery easy. Allows user to refold for transport or storage.
- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Roomy 19" between arms is easy to use.
- Wide base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Steel frame is rust resistant and chrome plated finish offers greater durability.
- Includes pail, lid and splash guard
- Height adjusts 16.5" - 22.5"
- 350 lb weight capacity.
- Lifetime limited warranty.

### Models:

- G30213-4F (Folding Frame - Painted)
- G30213-4 (Painted)
- G30211-4 (Chrome Plated)



Model G30210-4F

## 3 IN 1 ALUMINUM COMMODE

---

- Aluminum frame is lightweight and easy to transport.
- Folding frame design assembles in seconds without tools and reduces storage space by 35% while making delivery easy. Allows user to refold for transport or storage.
- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Full 19" between arms is easy to use.
- Wide base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Includes pail, lid and splash guard
- Height adjusts 17.5" - 21.5"
- 250 lb. weight capacity.
- Lifetime limited warranty.

### Model:

- G30210-4F (Folding Frame - Aluminum)



Model G30214-2

## HEAVY DUTY 3 IN 1 COMMODE

---

- Extra wide frame is easy to use with a full 21.12" between arms.
- Large base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Height adjusts 18" - 22"
- Weight capacity of 400 lbs.
- Lifetime limited warranty.

### Model: G30214-2