



Ace Medical Supply, LLC.

Better Equipment! Better Life!


www.acehomemedical.com

Short-Wave Wheelchair Cushion



Simple Secure Attachment

Two Elastic hook and loop straps located at the back and seat help anchor the cushion to the chair.



Unique segmented fins are specifically engineered to address the effects of shearing that are typical of institutional seating.

Patent-pending design has two distinct zoned orientations to reduce shear in multiple directions for specific seating areas.

(Use of cushion also helps to reduce pressure.)



Fluid-repellant Cover

Removable zippered cover features non-slip top fabric and durable non-skid bottom. Wipes clean with neutral suds or approved cleaning agents.

Item	Description	QTY
SWAVE - 04	Short-Wave Cushion W/ cover	4/cs
SWAVE - 01	Short - Wave Cushion W/Cover	1/cs

Warranty: 18 Months
 Weight Limit: 300 lbs
 Cushion Width: 18"

Designed to fit 16" or 18" deep seat

Ace Medical Supply, LLC. • 7515 N. Milwaukee • Niles, IL 60714

PH: 847-647-8590 • Fax: 847-647-7808