

Guardian® Bath Safety

Guardian® offers a complete line of bath safety products to help meet all of your customers' needs. With a commitment to delivering the highest level of quality, each Guardian product is engineered for safety and durability. These bath safety products offer the ideal choice for comfort and convenience.

GUARDIAN EASY CARE® SHOWER CHAIR/STOOL

A lightweight and economical choice for those who require a basic chair with height adjustability and an extra wide seating surface.

- Sturdy, lightweight, anodized aluminum frame resists rusting.
- Legs adjust in 1" increments with just a push of a button to ideal height between 13" and 20".
- Available with a back for extra support and leverage, or backless for easier access in bathing.
- Textured, blow-molded polyethylene wide track seat and backrest with built-in hand-holds provide extra safety.
- Rubber feet grip even wet surfaces.
- Available in attractive retail packaging.
- Weight capacity 250 lbs.

Model:

- G30400-4 (Chair assembled-Gray)
- G30401-4 (Stool assembled-Gray)
- G30402 (Chair unassembled-Retail-White)
- G30403 (Stool unassembled-Retail-White)



Model G30402



Model G0403



Model G30406B

GUARDIAN *Plus+* BARIATRIC SHOWER CHAIR

- Frame is reinforced to support weights up to 500 lbs.
- Extra wide seating surface accommodates larger users.
- Seat height adjusts from 15"-19" for a comfortable fit.
- Rubber tips grip bathtub surface.

Model:

- G30406B



Model G98027

GUARDIAN *Select* PADDED SHOWER CHAIR

- Extra thick padding on seat and backrest provides additional comfort for long-term use.
- Extra deep 16 1/2" seat for greater stability and comfort.
- Legs adjust in 1" increments to ideal height between 18" and 22".
- Rubber tips grip bathtub surface.
- 300 lb. weight capacity.

Model:

- G98027