



Model G98308A

NON-PADDED TRANSFER BENCH

- Textured, molded backrest and seating surface offer a non-slip surface with holes for water drainage.
- Lightweight, anodized aluminum frame is height adjustable in 1/2" increments.
- Backrest attaches for right or left facing tubs easily without tools.
- Side arm and seat back handle provide extra stability and leverage when transferring.
- Weight capacity 300 lbs.

Model:
G98308A (Assembled)



Model G98338F

PADDED TRANSFER BENCH

- Vinyl covered foam backrest and seating surface offer a highly resilient cushioned surface for comfort.
- Lightweight, folding aluminum frame is height adjustable in 1/2" increments.
- Suction cup footpieces secure the bench to the tub floor.
- Backrest adjusts for right or left facing tubs.
- Side arm and seat handle provide extra stability and leverage when transferring.
- A cost effective solution with the added feature of padding.
- Weight capacity 300 lbs.

Model:
G98338F (Assembled)



Model G98309

GUARDIAN^{Select} NON-PADDED TRANSFER BENCH

- Suction cup footpieces secure the bench to tub floor
- Lightweight, anodized aluminum frame is height adjustable in 1/2" increments to accommodate most tubs.
- Backrest adjusts without tools for right or left facing tubs.
- Built-in soap dish and hand-held shower slot for user convenience.
- Weight capacity 300 lbs.

Model:
G98309 (Assembled)
G98309G (Unassembled-Retail)



Model G98013A

PADDED TRANSFER BENCH WITH COMMODE OPENING

- A commode opening affords the user greater comfort and better access for cleaning.
- Thick, high density foam cushions and backrest for comfortable seating that won't bottom out.
- Suction cup footpieces secure the bench to tub floor.
- Lightweight, folding aluminum frame is height adjustable in 1/2" increments to accommodate most tubs.
- Weight capacity 300 lbs.

Model:
G98013A (Assembled)



Model G98310B

GUARDIAN^{Plus+} BARIATRIC TRANSFER BENCH

- 500 lb. weight capacity accommodates larger users.
- Textured, molded backrest and seating surface offer a non-slip surface with holes for water drainage.
- Aluminum frame is height adjustable in 1/2" increments.
- Backrest adjusts without tools for right or left facing tubs.
- Side arm and seat handle provide extra stability and leverage when transferring.
- Safely provides assistance to those who have difficulty stepping over bathtub walls.
- Backrest and armrest installed with detent buttons for tool free assembly.

Model:
G98310B (Assembled)